

LOC441150-similar to RIKEN cDNA 2310039H08 Gene

Reference: PTXBC060325

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to RIKEN cDNA 2310039H08
GENE ID	441150
ORIGIN SPECIES	Human
ACCESSIONS	BC060325
VECTOR	pDONR223
NCBI SYMBOL	LOC441150
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcgtcccgagtcgccaatgctcgccccggcctctgctcagcttgggttacctggcgagctcctgagctggtccagcagggccaggaaactccggcctggagaaacgccacatcgccgatccacggcggaaccacagcgtcccggtgccgaggagccca
agccctgggagggcgccgctttggctgagtccttccccctccgacctcaggatagggaacggccccggcgagcctggcttgggtgaggcagcgactgcgcttctcatggcatacagtgggccccctga