

MCOLN3-mucolipin 3 Gene

Reference: PTXBC060765

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TRP-ML3; TRPML3; mucolipin-3; mucolipin 3		

For research use only.

GENE DESCRIPTION	mucolipin 3
GENE ID	55283
ORIGIN SPECIES	Human
ACCESSIONS	BC060765
VECTOR	pDONR223
NCBI SYMBOL	MCOLN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagatcctgagtgatgtgtagtagctgagctctcatgaagaggaaaatcgctgcaattttaaccagcaaacatctccatctgaggagcttctattagaagaccagatgaggcgaaaactcaaatTTTTTcatgaatccctgtgagaagtctgggctcgaggtagaaaaccatggaa
acttgcatacaaaattcaaaaatgcaatggtagactatccagctggcttatttgggctaagtaaccagatgggtgtagctttcaagggaagagaatactatagcattcaaacaccttttccataaaggatataatggaccgaatggatgacacatatgcagtgtacacacaaatgacgtgtatgat
cagttaatcttcgagtaaacagctacttgcagctatacaatgtctcgttgggaatcatgcttatgagaacaaagggtaccaaagcaatctgctatggcaatctgtagcagcattatacaagcgaggaaacatctaccctggaaatgatactttgacatcgatccagaaatgaaactgagtggttctt
tgtggagccagatgaaccttttcaattgggacaccagcagaaaataaactgaacttaacactggacttccagactcctaacagtgaggcttcagtttaaaactgaaagccattatctgcagacagttcgtcatcaagaacctccctgactgttatgactttactctgactataacatttgacaaca
aggcccatagtggaagaattaaaataagtttagataatgacattccatcagagaatgtaaagactggcatgtaictggatcaattcagaagaacactattacatgatgatctttgatgcctttgtcattctgacttgccttgcatttaatacctctgcattagatctgtgattagagacttcagctca
gcaggtagggaacgttgccttctag