

TMEM154-transmembrane protein 154 Gene

Reference: PTXBC060775

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 154		
For research use only.			

GENE DESCRIPTION	transmembrane protein 154
GENE ID	201799
ORIGIN SPECIES	Human
ACCESSIONS	BC060775
VECTOR	pDONR223
NCBI SYMBOL	TMEM154
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggctccccgcgcagccctagctctgcctggtgatcgctcgtcccgctcgccggggtaattatgaggaattagaaaactcaggagatacaactgtggaatctgaaagaccaaataaagtactattccaagcacatttgcagtgaccatcaagaacattaaatgcaaat
ataaattctaccaactttgctccgatgaaatcagttagagtttatactgatggtttaatccattgattttattgctccttactttttattcgtggtattccttgcacatactataaaagaaaaagaactaaacaagaacctctagccaaggatctcagagtgctttacagacatatgaactggga
agtgaacacgtgaaagtcctatttttgaggagatacaccctctgttatggaattgaaatggaagagcttgataatggatgaacagcatgaatagaatgccgacttgaatgtttacctaccttgaggagaaggaatcaatcacaaccaagtgaacgtgaatcctaa