



Reference: PTXBC060810

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPR172B; GPCR42; PAR2; RBFVD; RFT1; RFVT1; hRFT1; solute carrier family 52, riboflavin transporter, member 1; G protein-coupled receptor 172B; G-protein coupled receptor GPCR42; PERV-A receptor 2; porcine endogenous retrovirus A receptor 2; solute carrier		

GENE DESCRIPTION	G protein-coupled receptor 172B
GENE ID	55065
ORIGIN SPECIES	Human
ACCESSIONS	BC060810
VECTOR	pDONR223
NCBI SYMBOL	GPR172B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]