



Reference: PTXBC060820

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZBP-99; ZNP-99; zinc finger protein 281; GC-box-binding zinc finger protein 1; ZNP-99 transcription factor; transcription factor ZBP-99; zinc finger DNA-binding protein 99		

For research use only.

GENE DESCRIPTION	zinc finger protein 281
GENE ID	23528
ORIGIN SPECIES	Human
ACCESSIONS	BC060820
VECTOR	pDONR223
NCBI SYMBOL	ZNF281
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]