

PTPLB-protein tyrosine phosphatase-like (proline instead of catalytic arginine), member b Gene

Reference: PTXBC060839

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PTPLB; very-long-chain (3R)-3-hydroxyacyl-CoA dehydratase 2; protein tyrosine phosphatase-like (proline instead of catalytic arginine), member b; protein-tyrosine phosphatase-like member B; very-long-chain (3R)-3-hydroxyacyl-[acyl-carrier protein] dehydra		

For research use only.

GENE DESCRIPTION	protein tyrosine phosphatase-like (proline instead of catalytic arginine), member b
GENE ID	201562
ORIGIN SPECIES	Human
ACCESSIONS	BC060839
VECTOR	pDONR223
NCBI SYMBOL	PTPLB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcagtgccggcgactgcagcagcgaagggaatggggcgccggggcaggccggggccggggcaggccagcgccagcgggaagaagaaggccggggccctggccacggcgctacgtgtcatctacaatgtgtgatgcagccgggtggctgttatagcggttggtctgtccgagcatacctggcctaaaggtagctacacatagcctttattatcaattgaaaagccttgaattctttcaaaactggagccttattggagattttacattgtgctataggaattgtccatctctgttgctgactcttccagggtgatgtcaagagttttctaataatggcagtaacaatagcgtcaaaagggtacagagtgaaacagagtgctcctgtttgttattgcatggagatcacggaaatcatcgttactccttttatacatcagctctataaaccatctgccttacctcatcaaatggggcagggtacacactttcattgtgctgtacccaatgggagtgctcaggagaactgctcaaatatgcagctctgccctgttcagacaagctggcctatatccatcagtttacccaacaatacaattctctttgactactatgcattctgattctataatgatctcctacattccaattttcccccagttatacttccacatgataccagagagaagaagatcctttctcactgaagaacaagaatttgatag