

SOCS3-suppressor of cytokine signaling 3 Gene

Reference: PTXBC060858

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATOD4; CIS3; Cish3; SOCS-3; SSI-3; SSI3; suppressor of cytokine signaling 3; STAT-induced STAT inhibitor 3; cytokine-inducible SH2 protein 3		

For research use only.

GENE DESCRIPTION	suppressor of cytokine signaling 3
GENE ID	9021
ORIGIN SPECIES	Human
ACCESSIONS	BC060858
VECTOR	pDONR223
NCBI SYMBOL	SOCS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcaccacagcaagttcccgccgggatgagccgccccctggacaccagcctgcctcaagacctcagctccaagagcgagtaccagctggtgtaacgcagtgccgaagctcagagagcggcttctactggagcgcagtgaccggcggcgaggcgaacctgt
gctcagtgccgagcccggcacttctgatccgacagctcgaccagcgccacttctcagctcagcgtcaagaccagctctggaccaagaacctgcgcctcagtggtgagggggcagcttctctcagagcagatccccggagcagcagcccgtgccccgttcgac
tgcgtctcaagctggtgcaccactacatgccgccccctggagccccctctccctcgcacctactgaacctctccgaggtgcccgagcagcgtctgccagccactcctggagtcccccagaagagcctattacatctactcggggcgagaagatccccctggtgtg
agccggccccctctccaagctggcactctcagcatctctgctggaagacgtcaacggccacctggactctatgagaaagtcaccagctgccggggccctcgggagttcctggaccagtagcctccctctaa