



Reference: PTXBC060861

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf39; AGO61; GTDC2; MDDGA8; protein O-linked-mannose beta-1,4-N-acetylglucosaminyltransferase 2; glycosyltransferase-like domain containing 2; glycosyltransferase-like domain-containing protein 2; protein O-linked mannose N-acetylglucosaminyltransferas		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 39
GENE ID	84892
ORIGIN SPECIES	Human
ACCESSIONS	BC060861
VECTOR	pDONR223
NCBI SYMBOL	C3orf39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]