

NDUFA3-NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 3, 9kDa Gene

Reference: PTXBC061644

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CI-B9; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 3; NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 3, 9kDa; NADH-ubiquinone oxidoreductase B9 subunit; complex I B9 subunit; complex I-B9; NADH:ubiquinone oxidoreductase subunit A3		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 3, 9kDa
GENE ID	4696
ORIGIN SPECIES	Human
ACCESSIONS	BC061644
VECTOR	pDONR223
NCBI SYMBOL	NDUFA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcgagagtcggcgcttctcaagaatgcctgggacaggagccagtgctgctgtgctctctgctgctgggggctcgtgttaattctgccccattgagccctacttcaagtactccgtcatgaacaaggccacgccctacaactaccagtgcccgctcgatgatgg
gaacatgcccagctgcccagccaccccaaggacctcaggggcccgagctgagtgctgaagaactgtga