

C3orf62-chromosome 3 open reading frame 62 Gene

Reference: PTXBC061697

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C3orf62; chromosome 3 open reading frame 62		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 62
GENE ID	375341
ORIGIN SPECIES	Human
ACCESSIONS	BC061697
VECTOR	pDONR223
NCBI SYMBOL	C3orf62
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgcattacataaagacatgctcgttttaggagaaatgctgaaaaactaagaagatgcagaaaggaactgactgcagccattgaccggccttgaaggagtagttattcccagagtgccacagccagcagagcctggaactgagcggcggccgctctctctcgtgccgtgc
acaggctcctctgcagaaagacatcctctgagcctgctctctgctcctcttctgctgctccatgtgctcctgagaatgagaacctgcttgcacaacacatgccccgtaaatgcaaaaccacatgctctgccccgagagaaaaccttaaccagcaaggaaaatgattgatg
cattcctccatttggcaccctgaaagagagctctggagaactgcaggagaggggaaaaactggagaaaaaataaggaagatggagagagattgaaaggctgactcaaacatgccactcaacaatttagccaagaggtcacaaggatctgctgatgatgaccatacaag
catccgaactatgaaagattgctggaataatgaaatgaaaacgaattgaaccacatgtgtgctcattccaagattcaccctcaagaggaaagcctggcctgctcatgacaagagccctcagaaggccacagatgctgacctggcagcctcaaacaggctttgatgataaa
tattgtgagactgtctgactggaagagactacaatgtgatgactctttaaataccaatgagtaa
```