

SYPL1-synaptophysin-like 1 Gene

Reference: PTXBC061887

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H-SP1; SYPL; synaptophysin-like protein 1; pantophysin; synaptophysin like 1		
For research use only.			

GENE DESCRIPTION	synaptophysin-like 1
GENE ID	6856
ORIGIN SPECIES	Human
ACCESSIONS	BC061887
VECTOR	pDONR223
NCBI SYMBOL	SYPL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcccaacatctacttggctgccagcggatcagtcgactcggcagaggatgtccggctccagatcaacctcaaccgctcaaggagccactcggctcctcaaggctcctcagtgagggttttaaggccaaacagaaattcaagt
gaattgtcctcctgcagttactgagaataaaactgttacagctacttttggatcattcagggtgaatgagggatcattcagccacctccaggtgtaaacatattgtagttaaattggaagattacgtcctcataggcgattactcttctgcacaattctatgttacctttgcagctttgtgttcct
gtactgcattgctgccctctgctttatgttggtctacacagagctgtatctggatagctgtaaacctcctatgatactttgtgttacactttgtgccactttttgtggttggtgagcacttcagcctgggctaaagctctgacagatattaaaatgctactggtcacaatatattgatgaactccgc
ctgtgaagaagaagcagtagctgttactttggctctgtgaccagatggggtccctaaatgtatctgtgatatttggctttctaaatgatactctggggaggaatgcttgggtgtgtacaaggagaccagcctacacagtcctcaatacatctgccctcatagccaaggaggtattcca
cctcctaccggaataaa