

C14orf49-chromosome 14 open reading frame 49 Gene

Reference: PTXBC061899

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf49; C14orf139; LINC00341; NCRNA00341; NET53; Nesp3; nesprin-3; long intergenic non-protein coding RNA 341; nuclear envelope spectrin repeat protein 3; spectrin repeat containing nuclear envelope family member 3		

For research use only.

GENE DESCRIPTION	chromosome 14 open reading frame 49
GENE ID	161176
ORIGIN SPECIES	Human
ACCESSIONS	BC061899
VECTOR	pDONR223
NCBI SYMBOL	C14orf49
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggaaacgctcttaattgtgcacactgtcctgtgtctgcagagcctgtcgtcagcacaagactttggagcagcttttgagccctgcagaggaaagctcttgacctccaggtcagggccaagccgagaaggggttcagcgggaccttctggaaaacaggcccagctctcaaggtgcaggtgaggggtctgtga