

OBSL1-obscurin-like 1 Gene

Reference: PTXBC061909

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	obscurin-like protein 1; obscurin like 1		
For research use only.			

GENE DESCRIPTION	obscurin-like 1
GENE ID	23363
ORIGIN SPECIES	Human
ACCESSIONS	BC061909
VECTOR	pDONR223
NCBI SYMBOL	OBSL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggggccccatcgccgctggtgctgcccgcacccagccctcagacggggcgagtttcagtcgctgctggagatgagtgctgctacttcactgtcaccatcacagatgtctcctcgatggtatgccagcggcaaggtgtatgtggcagcagtcgcctggagcgtgtgg
tgctga