

HIST1H2AA-histone cluster 1, H2aa Gene

Reference: PTXBC062211

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2AA; H2AFR; TH2A; bA317E16.2; histone H2A type 1-A; H2A histone family, member R; histone 1, H2aa; histone H2A/r; histone cluster 1, H2aa; histone cluster 1 H2A family member a		

For research use only.

GENE DESCRIPTION	histone cluster 1, H2aa
GENE ID	221613
ORIGIN SPECIES	Human
ACCESSIONS	BC062211
VECTOR	pDONR223
NCBI SYMBOL	HIST1H2AA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggacgagggaagcaggaggaagaaagcagcgccaagtctaagtcgtcttctagagcgggtttgcagtttccgtaggcggatccatcgtctcttctgaagggaactatgcagagcggatagggcaggcgcaccagtgtattggcgcagtgtagtatctcaca
gcagaaatccttgagctggcaggcaatgcgtctcgataacaaaaaactcgcaattatccccgccactgcagctagcgaatgatgaggaactcaataagctttggcgccgtgaccaatgccaggcgaggctctgctaactcaggcagtgctgctgccaagaag
actgagagtcaccaccataaagcccaagcaagtaa