

DEFB119-defensin, beta 119 Gene

Reference: PTXBC062212

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DEFB-19; DEFB-20; DEFB120; DEFB20; ESC42-RELA; ESC42-RELB; beta-defensin 119; defensin, beta 120; defensin, beta 19; defensin, beta 20; defensin beta 119		

For research use only.

GENE DESCRIPTION	defensin, beta 119
GENE ID	245932
ORIGIN SPECIES	Human
ACCESSIONS	BC062212
VECTOR	pDONR223
NCBI SYMBOL	DEFB119
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaactctttacctgtttcttgcctccttctggccatagaagaaccagtgtatcaggcgaacgccatccttcgatgcattgggtaacagtggaattgtaggccctcttgcaaaaagaacgaacagccctacctctattgcagaaattgtcagtcctgctgctccagtcctacatgagga
taagcatttctggcaagaggaaaataccgactggctcttatgagaagcagtgccaagactacctga