

C4orf7-chromosome 4 open reading frame 7 Gene

Reference: PTXBC062213

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C4orf7; FDC-SP; follicular dendritic cell secreted peptide; FDC secreted protein; follicular dendritic cell secreted protein		

For research use only.

GENE DESCRIPTION	chromosome 4 open reading frame 7
GENE ID	260436
ORIGIN SPECIES	Human
ACCESSIONS	BC062213
VECTOR	pDONR223
NCBI SYMBOL	C4orf7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaaagttctctctgatcacagccatctggcagtggtgtgtgttccagtcctcaagaccaggaaacgagaaaaagaagatcagtgacagcgatgaattagcttcagggtttttgtgtccctaccatataccattcgcccaactccaccaattccattccaagattccatggttagacgtaattttctattccaatacctgaatctgccccctacaactccccctctagcgaaaaagtaa