

LACRT-lacritin Gene

Reference: PTXBC062217

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	extracellular glycoprotein lacritin		
For research use only.			

GENE DESCRIPTION	lacritin
GENE ID	90070
ORIGIN SPECIES	Human
ACCESSIONS	BC062217
VECTOR	pDONR223
NCBI SYMBOL	LACRT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaattaccactctctcttggcagctgtacagggggccctgggtatgtgaagatgcctctctgactcgacgggtgctgatctgccagggaagctgggaccttaagcctaataagagatctcaggtccagcagaaccagcttcaccccagagacaaccacaacagccca
ggagacttcggcgagcaggttcaggggacagccaaggtcacctcaagcagggcaggaactaaacccctgaaatccatagtgagaaaagtatcttactaacagacaagcccttgcaaaagcaggaaaaggaaatgcacggaggcgtgccaggtggaaaacaattcatcgaaaatg
gaagtgaattgcacaaaaattactgaagaaatcagtcatttaaaacatgggcatga