

SCGB2A1-secretoglobin, family 2A, member 1 Gene

Reference: PTXBC062218

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LPNC; LPNC; MGB2; UGB3; mammaglobin-B; lacryglobin; lipophilin C; mammaglobin-2; secretoglobin family 2A member 1		

For research use only.

GENE DESCRIPTION	secretoglobin, family 2A, member 1
GENE ID	4246
ORIGIN SPECIES	Human
ACCESSIONS	BC062218
VECTOR	pDONR223
NCBI SYMBOL	SCGB2A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagctgctgatggtcctcatgctggcggccctcctcctgactgctatgcagattcgtgctgcaaacctcctggaggacatgggtgaaaagaccatcaattccgacatatctatacctgaatacaaaagagcttctcaagagttcatagacagtgatgccgctgcagaggctatgggaaatt
caagcagtggttctcctcaaccagtcacatagaactctgaaaaacttggactgatgatgcatacagtgtagacagcatttgggtgaatatgaagagtaattaa