

FNDC5-fibronectin type III domain containing 5 Gene

Reference: PTXBC062297

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FRCP2; irisin; fibronectin type III domain-containing protein 5; fibronectin type III repeat-containing protein 2; fibronectin type III domain containing 5		

For research use only.

GENE DESCRIPTION	fibronectin type III domain containing 5
GENE ID	252995
ORIGIN SPECIES	Human
ACCESSIONS	BC062297
VECTOR	pDONR223
NCBI SYMBOL	FNDC5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgtcctcatcaggaggtgaacaccaccaccgctcatgtgccctctgggacctggaggaggatagcaggtacatagtcacgtgcaggccatctccattcaggccagagcccagccagcagcctgtgctctcaagaccccgctgaggctgagaagatggcctcaaga
acaaagatgaggttaacctgaaagagatggggaggaaccaacagctgcggacagcggaggtgctgatcatgctgtgctctgttcattgcccctctctgcccagctatgacatcatcaaggacaatgaacccaataacaagaagaaaaaccaagatgca
tcagaaccagcacaccagagcaccagggcgggggcttccgcagcaaggtgagggcaagacctgggcctgggtggccaccctgtgcctcatgctctggtaa