



Reference: PTXBC062320

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EXLD1; HsT19883; hyaluronan and proteoglycan link protein 3; extracellular link domain containing, 1		

For research use only.

GENE DESCRIPTION	hyaluronan and proteoglycan link protein 3
GENE ID	145864
ORIGIN SPECIES	Human
ACCESSIONS	BC062320
VECTOR	pDONR223
NCBI SYMBOL	HAPLN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1