



Reference: PTXBC062323

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf25; C21orf258; TMEM24L; C2 domain-containing protein 2; TMEM24 related gene; C2 calcium dependent domain containing 2		

For research use only.

GENE DESCRIPTION	C2 calcium-dependent domain containing 2
GENE ID	25966
ORIGIN SPECIES	Human
ACCESSIONS	BC062323
VECTOR	pDONR223
NCBI SYMBOL	C2CD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]