



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

IL27-interleukin 27 Gene

Reference: PTXBC062422

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL27-A; IL-27; IL-27A; IL27A; IL27p28; IL30; p28; interleukin-27 subunit alpha; IL-27 p28 subunit; IL-27 subunit alpha; IL-27-A; interleukin-30; interleukin 27		

For research use only.

GENE DESCRIPTION	interleukin 27
GENE ID	246778
ORIGIN SPECIES	Human
ACCESSIONS	BC062422
VECTOR	pDONR223
NCBI SYMBOL	IL27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagcagcagcagcagcacttggctggcgccgcagcgcgtgtgtctgtcttccctgtcctctggttcaacgtgctgtctggaggttcccaaggccccaggaggagccccagctgagcctgcagcagctgcggaggaggttcacagctacgctgtcatctccgccaagctgctcgcgcg
agggttccggggccaggccccaccgttttgcggaatctcactgtccaggagtgtaactctgtacctctgccccctgggagagcagctcctgatgtttccctgaccttccaggcctggcgccgcctctgtacccggagcgtctctgttcatctccaccacgtttacgccccttcatgcccctgtggtg
gagggcctggggaccacaggcccgctgaccacaacttgagagagtgagcagctgtggggccatgagggttgagactccgcgatctgcagcggcaccctccgcttccaggtgtctgctgcacgattcaacctcccggagagaggaggagaggaagagagaggaggagaggaggagaa
ggggcctgtccacagcagctgggcagcgccttacagggcccgcccgagtgctctgccccagctccttccacatccgcctgtgctactctgtggagctgctgtatctatctgggcccgtgcggggagttgtctgctgtctgcacaaggctgggcactacgtgtggcccttgggttcccaac
attgaacccccaacccctaa