

CCDC150-coiled-coil domain containing 150 Gene

Reference: PTXBC062449

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	coiled-coil domain-containing protein 150; coiled-coil domain containing 150		
For research use only.			

GENE DESCRIPTION	coiled-coil domain containing 150
GENE ID	284992
ORIGIN SPECIES	Human
ACCESSIONS	BC062449
VECTOR	pDONR223
NCBI SYMBOL	CCDC150
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgagaatgtgtggtctcagatgaatgcaaggtgctctggagaaagtacaatagacttggcggaggattcagagattgcaggcctcaagaaagaaaggatcacaacagaggggtgcagaagctggaagctgaagtggaccagtggcaggccaggatgctgt
catggaggaccagcacaacagtgagattgaatctctacaaaagcttaggtgtagctagagaagacaacaggaaacttgctatgagcttggaacaagctctccagacaaataatcatctgcaacaaagctagatcacattcaagcaattggaaagcaagaacttgagcgacagaat
ttggaaccttcaagaccggatgactgaagagtcacaagtggagcagaattgcatgctgaacgcatagaagctctagaagaagcagtttcaaacggagagagaaactacaagaaagtgccacaacgggaagtgccgagctgaagaaagcccttgatgaagcctaactcagatcagt
ggaagtgtcccgaccacccgagagctgcgacagaaactgcagagctagaaaaatactagaagaagtaacaaggagaaaaataagaatcaaaagacccaattaagctccactgtcagctaaggcgaataatgctcagaatataaggaatgaaggtgtatgggaaacctctctca
cttctcgataacctgtga