

FLJ40125-hypothetical protein FLJ40125 Gene

Reference: PTXBC062452

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein phosphatase, Mg ²⁺ /Mn ²⁺ dependent 1N (putative)		
For research use only.			

GENE DESCRIPTION	hypothetical protein FLJ40125
GENE ID	147699
ORIGIN SPECIES	Human
ACCESSIONS	BC062452
VECTOR	pDONR223
NCBI SYMBOL	FLJ40125
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacctgcatcctgtctgttccctggggccctaggcctctgaggaggcgtcaggaggagctagcactggacgcagccctgggctgcagaatcgtgaactgtgtgcctctgctcagaagccccagcctgaacacagtttccaggactctggcctcagaggacatccagatt
tacctctggggaggggctggactgcaaggccactgtcattgtgaagtgttctcagatctgccaggctcagaagagtgccggagagaagggcaggatgggctgggaagtccaacccacgcatttgggctcagcctggacatggaggcctga