

C5orf33-chromosome 5 open reading frame 33 Gene

Reference: PTXBC062567

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C5orf33; DECRD; MNADK; NADKD1; NAD kinase 2, mitochondrial; NAD kinase domain-containing protein 1, mitochondrial		

For research use only.

GENE DESCRIPTION	chromosome 5 open reading frame 33
GENE ID	133686
ORIGIN SPECIES	Human
ACCESSIONS	BC062567
VECTOR	pDONR223
NCBI SYMBOL	C5orf33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgctggcagcgagtaagcttggacagacttaaccagttatagggttaaacactgatccagaacggtctgagggtcatttatgctgcccgttcgatatcacattccttcagaagccttacagaagttctatcggtgagttcaggtggttgaggcagagaatcaggttatac
cttgaagggactggcataaaccctgtacctgtggaccttcacgagcagcagctaaagctgaatcagcacaatagagcccttaacattgaaagagctcatgatgaaaggtctgaggcttcaggaccccaactctgagtgagagcactaaatgaagttctattgggagagctgtcatcc
agggttcctactatgagattcagttgatgatggtccatgggaaaaacagaaggttcagggtcctcaatttgtgactggaacaggatcaaaaggcctgggtcattcaatattaacaggggtgcaactcaggtctgtagaagatgttttaattgcaaacgacaaggaatttgagttcctcattga
acagagaattggtagagaaagtaacaatgaatataatgaatcactgctctacagtcgggaagaacaaaatactttcagttatcagaaaccaatagcaaatagagtttctcaagcagtcgagcgtgtgttctctcaaggtttgttcgttctcgttgggatgcctgtatggttgg
atggaggaactcttttgagtttaatgatggtgcaattgcttcgatgatcaataaagaagatgagcttgaactgtgcttctgaacagtga