

C1orf104-chromosome 1 open reading frame 104 Gene

Reference: PTXBC062571

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf104; RUSC1 antisense RNA 1 (non-protein coding); RUSC1 antisense gene protein 1; RUSC1 antisense RNA 1		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 104
GENE ID	284618
ORIGIN SPECIES	Human
ACCESSIONS	BC062571
VECTOR	pDONR223
NCBI SYMBOL	C1orf104
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagccgggaggggtcagagaacgcggctgcgtctggaatcctgaagggggcgggggcctggaagaggtccggggccagaatggacctcaggtccctgctgccgaatctggccctgccctccagctaccccttactccagaggaagggaaccaggagaccaccca
gatgccctgaaaggaggtggaggtgggggtggggcaatagcagagttgtctggtgaatgcaggaagggtggggcctggagaggaaggaatggagcagccgtcagctgtgcaactccgcatctgctggctgctcggcggtccagcctgcgcctcggccctcggaaac
agcgggtatgtccattctctcgcattcgtctccatcgttcagccacacgcaccccttctctattcattctctcgtcctcctgtcctcggcccccctcgtgcccaggctcctggtctacacccatggggcgagccctcctcaccagggtcctcctggaaccgctcggccttgggcttgcctg
cggctccacggtcccccggcgggggcacagagcgggaggggggacactagcacaacctaccctgcgtgcgcagctgcacccctgcgcgatggcggtggcgtagctcagaccgccccagcgttttagcgtctttgtacccacctagagggtttgatatacctaa