

TMEM71-transmembrane protein 71 Gene

Reference: PTXBC062592

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 71		

For research use only.

GENE DESCRIPTION	transmembrane protein 71
GENE ID	137835
ORIGIN SPECIES	Human
ACCESSIONS	BC062592
VECTOR	pDONR223
NCBI SYMBOL	TMEM71
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtaccgaatctcaactgatgtcaacaccagtagcaagttctccaggttgaaagagaatatgctggagagctgtctccacgtgcatttcccaagtttcacgtgtattccctggatgggtaccattctttgaatgcggctccatagatccctgacaggtccactatactgtcgcg
aagtccagactcccaaatggctactatattggactgaagacagcttctgtgcgacaaagatggcaacataactctgaacccatccagaccagcgttatgtataaggagaacttagttagtaccccaatcttggctgcagtggaagtatcttgggtgacatcaactcttccaagtga
gacaactggttgaaagggaccagaggttgacacagaccattgcaatggaatgcagatgatttagactgttctctctgactgatgactggagtcaggggaagatgaatgcagagtctgtgatcacctctctccagccacatcatctcagcctcctggaggaaactcccatagctgt
ctcttcagtcctcagacgttctgaacgtttccaagagaatagttcggatcagaaaccaggtgtgtgcaagaggtcttcttcaggcaatcctgctgtgtgcttaactattctgcatgtgcaagatgggttatgggagaataatagccagtgcttcacatgctcattgatgataact
gtagcttatgtgaatcattgttctcagccttgccagctatttcaaaaccactgcctgtgctggttgcaaaattga