

C22orf39-chromosome 22 open reading frame 39 Gene

Reference: PTXBC062599

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0545 protein C22orf39; chromosome 22 open reading frame 39		
For research use only.			

GENE DESCRIPTION	chromosome 22 open reading frame 39
GENE ID	128977
ORIGIN SPECIES	Human
ACCESSIONS	BC062599
VECTOR	pDONR223
NCBI SYMBOL	C22orf39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggacggcagcggcctggcagccgcccctgcgagggcctaccgcccagtggaagctctgccgcagcggcaggcacttctacaccactactacgtccacggcgagcggcggcctgcgaacagtggcagcggcagctggccagctgccgcgactgggaggagc
gccggaacggcggagccagcaatccctctgtgagagcgagcgggcacgagtcgggctgcacggaagcacatcctggtgtgggccccgaggcagacccccctccagactggcatctccctctgccacaggagaaggacgagtgta