

PKN2-protein kinase N2 Gene

Reference: PTXBC062620

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PAK2; PRK2; PRKCL2; PRO2042; Pak-2; STK7; serine/threonine-protein kinase N2; PKN gamma; cardiolipin-activated protein kinase Pak2; protein kinase C-like 2; protein-kinase C-related kinase 2; protein kinase N2		

For research use only.

GENE DESCRIPTION	protein kinase N2
GENE ID	5586
ORIGIN SPECIES	Human
ACCESSIONS	BC062620
VECTOR	pDONR223
NCBI SYMBOL	PKN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaaaaagtaaagccaccatttatacctaccataagaggacgagaagatgtagtaattttgatgatgaatttacctcagaagcacctattctgactccactcgagaaccaaggatacttcggaaggagcaggaaatgttcagagattttgactacattgctgattggtgtaa