



Reference: PTXBC062625

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AEZ; A WMS2; ZIP4; zinc transporter ZIP4; solute carrier family 39 (zinc transporter), member 4; zrt- and Irt-like protein 4; solute carrier family 39 member 4		

For research use only.

GENE DESCRIPTION	solute carrier family 39 (zinc transporter), member 4
GENE ID	55630
ORIGIN SPECIES	Human
ACCESSIONS	BC062625
VECTOR	pDONR223
NCBI SYMBOL	SLC39A4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]