

LOC400590-hypothetical LOC400590 Gene

Reference: PTXBC062632

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC400590
GENE ID	400590
ORIGIN SPECIES	Human
ACCESSIONS	BC062632
VECTOR	pDONR223
NCBI SYMBOL	LOC400590
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcattttgtttatgagcaaggtgggtctcagagggtatcgccgatcagagggcgatgaagttctagatccattgagacaagctctagacagtagcatgcagtcgccacaactgtaccagcatcccgatctggcattccatgtttctgtctctgtggcctccacagtgcacaagcta
gcggtttacttggacctctacatctttcttttgcgttcagcctgcgattcgtttctctccacttggctctcatggtgcaggtttccaagaaatggcgctaaggccgagagttaa