



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ADAM33-ADAM metallopeptidase domain 33 Gene

Reference: PTXBC062663

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf153; DJ964F7.1; disintegrin and metalloproteinase domain-containing protein 33; ADAM 33; a disintegrin and metalloprotease 33; a disintegrin and metalloproteinase domain 33; disintegrin and reprolysin metalloproteinase family protein; ADAM metallope		

For research use only.

GENE DESCRIPTION	ADAM metallopeptidase domain 33
GENE ID	80332
ORIGIN SPECIES	Human
ACCESSIONS	BC062663
VECTOR	pDONR223
NCBI SYMBOL	ADAM33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttggagggcccgagagctcgggggaccccgctgctgctgctactactgctgctgctgctggccagtcgacggcgccgggggtgctcaaggacatataccttgggcagccagtcaccccgactgggtcctggatggacaacctggcgccaccgtcagcctggaggagccg
gtctcgaagccagacatggggcgctgctgcccctggagctgaagcccaaggactcctgctgagctggagaagaaccagtga