

WFDC11-WAP four-disulfide core domain 11 Gene

Reference: PTXBC062670

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein WFDC11; WAP11; protease inhibitor WAP11; WAP four-disulfide core domain 11		

For research use only.

GENE DESCRIPTION	WAP four-disulfide core domain 11
GENE ID	259239
ORIGIN SPECIES	Human
ACCESSIONS	BC062670
VECTOR	pDONR223
NCBI SYMBOL	WFDC11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcagcctcatgaagctctggataccatgctcatgacattctctgtacggtgctactgtctgtctgtggagaaatgaggagaaaagatatgacaggaagggaattgttactgaagaatgctgggaaagccaaatgtcaagaatgtaccaataagtggttctaaagcctttagatgtaa
agacaaaaattacacatgctgctggacctattgtggaacatctgctggataaacgtcgaaccagtgagattactaa