

RPL23-ribosomal protein L23 Gene

Reference: PTXBC062716

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L23; rpL17; 60S ribosomal protein L23; 60S ribosomal protein L17; ribosomal protein L17; ribosomal protein L23		

For research use only.

GENE DESCRIPTION	ribosomal protein L23
GENE ID	9349
ORIGIN SPECIES	Human
ACCESSIONS	BC062716
VECTOR	pDONR223
NCBI SYMBOL	RPL23
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgaagcgaggacgtggtgggtcctctgtgctgcaaaatccggatttcttgggtcttccggtaggagctgtaataattgtctgacacacaggagccaaaacctgtatcatctccgtgaaggggatcaaggagcgtgaacagactcccgtggtggtgacatggtgatggccacagtcaagaaaggcaaacagagctcagaaaaaggtacgtccagcagtggtcattcgacaacgaaagtcataccgtagaaaagatggcgtgtttctttatgtgaagataatgcaggagtcatagtgaacaataaaggcgagatgaaagggttctgccattacaggaccagtagcaaaaggagtgctgagactgtggcccccgaattgcatccaatgctggcagcattgcatga