

MGC72080-MGC72080 pseudogene Gene

Reference: PTXBC062733

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ASNSD; TS11; asparagine synthetase [glutamine-hydrolyzing]; TS11 cell cycle control protein; glutamine-dependent asparagine synthetase; asparagine synthetase (glutamine-hydrolyzing)		

For research use only.

GENE DESCRIPTION	MGC72080 pseudogene
GENE ID	389538
ORIGIN SPECIES	Human
ACCESSIONS	BC062733
VECTOR	pDONR223
NCBI SYMBOL	MGC72080
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcctggtcacggtgctgggaacctgctcatcctggccatcagccctgactccacccacacccccatgtacttctctccacactgtccttgctgacatcggtttcacctccaccacgggtcccaagatgattgtggacatcagctctcacagcagagtcctctatgca
ggctgcctgactcagatgtctctctttgccatttttggagcatggaagagagacatgctcctgagtgatggcctatga