

LCN6-lipocalin 6 Gene

Reference: PTXBC062746

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	epididymal-specific lipocalin LCN6; LCN5; UNQ643; hLcn5; epididymal-specific lipocalin-6; lipocalin 5; lipocalin 6		

For research use only.

GENE DESCRIPTION	lipocalin 6
GENE ID	158062
ORIGIN SPECIES	Human
ACCESSIONS	BC062746
VECTOR	pDONR223
NCBI SYMBOL	LCN6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcgctgctgctgctcttctgcttgggtctcgtgccaggccaggccgtgtgttgggaagactggaccctgagcagcttctggccctggtacgtgcttgcgggtggcctccgggaaaaggccttggcatggagaaggacatgaagaacgtcgtgggggtgtg
gtgaccctcactccagaaaacacctgcggacgtgctcctcagcacggcctggagggtgtgaccagagtgtcatggacctgataaagcgaactccggatgggtgtttgagaatccctcaataggcgtgctggagctctgggtgctggccaccaacttcagagactatgccatctc
tcactcagctggagttcggggacgagccctcaacaccgtggagctgtacagctgactgagacagccaggagccatggggctcttcaccaagtggagcaggagcctgggcttctcctgtcacagtag