

ABCA3-ATP-binding cassette, sub-family A (ABC1), member 3 Gene

Reference: PTXBC062779

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ABC-C; ABC3; EST111653; LBM180; SMDP3; ATP-binding cassette sub-family A member 3; ABC transporter 3; ABC-C transporter; ATP-binding cassette transporter 3; ATP-binding cassette, sub-family A (ABC1), member 3; ATP binding cassette subfamily A member 3		

For research use only.

GENE DESCRIPTION	ATP-binding cassette, sub-family A (ABC1), member 3
GENE ID	21
ORIGIN SPECIES	Human
ACCESSIONS	BC062779
VECTOR	pDONR223
NCBI SYMBOL	ABCA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtgctcaggcagctggcgtcctctctggaagaactacacctgcagaagcgaaggtctgtgacgtctggaactcttctgcattgctgtttctgggatcctcatctggctccgttgaagattcagtcgaaaaatgtcccaacgccaccatctacccgggcccagtc
atccaggagctgcctctgtttctaccttccctccgagagacacctgggagcttgctacatcccttctcacagtgcgctgccaagaccgtcactgagacagtgccgagggcacttggatcaacatgcgagtcgctggttccctccgagaaggacttggagactatagggtac
gacaactgctgtccagcgtgtggtggccgtgtgtctcgagcaccctcaaccacagcaaggagccctgcccgtggcgtgaaatcacctacggttcagttacacacggagaattacatgtggacccaacaggctccttttctgaaagagacagaaggctggcacactact
ccctttccgcttttccaaaccaggaaccaagggaacctacatccctgatggcggagaacctgtgtgagaagctcggttag