

SLC37A2-solute carrier family 37 (glycerol-3-phosphate transporter), member 2 Gene

Reference: PTXBC062990

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	glucose-6-phosphate exchanger SLC37A2; pp11662; solute carrier family 37 (glucose-6-phosphate transporter), member 2; solute carrier family 37 (glycerol-3-phosphate transporter), member 2; sugar phosphate exchanger 2; solute carrier family 37 member 2		

For research use only.

GENE DESCRIPTION	solute carrier family 37 (glycerol-3-phosphate transporter), member 2
GENE ID	219855
ORIGIN SPECIES	Human
ACCESSIONS	BC062990
VECTOR	pDONR223
NCBI SYMBOL	SLC37A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttctgtacaactacattggccaggacgggattgccagctccatagtgatgctgatcatctgtgggggctgtgcaatggccatacgcgctcaccactgtgtctctgtgacctggggactcacaagagcctgaagggaacgcaaaagccctgtccacggtcacggccatcat
tgacggcaccgggtccataggtgcggctctggggcctctgctgctgggtcctcctccacgggctggaacaatgtcttctacatgctcatctctgccgagctcctagcctgtgtgctcttggcgggttagtatacaagagatcttggcctggaaggtgtccctgagcagagcagcgg
gtataaagaatatga