

RBP7-retinol binding protein 7, cellular Gene

Reference: PTXBC063013

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRABP4; CRBP4; CRBP4V; retinoid-binding protein 7; cellular retinoic acid-binding protein 4; cellular retinoic acid-binding protein IV; cellular retinol binding protein 7; retinol binding protein 7, cellular; retinol binding protein 7		

For research use only.

GENE DESCRIPTION	retinol binding protein 7, cellular
GENE ID	116362
ORIGIN SPECIES	Human
ACCESSIONS	BC063013
VECTOR	pDONR223
NCBI SYMBOL	RBP7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgcgcacctcagcggtacttgaccctgctcagcagcgacaacttcgagggctacatgctggccctaggtattgactttgccactcgtaaaatagccaagtgtgctgaagccacagaaagtgattgagcagaatggggattctttaccatccacacgaacagcagcctaaggaacta
ctttgtgaaatttaaagttggagaagaatttgatgaagataacagagcgctggacaacagaaaatgcaagagtttggtatctgggacaatgacaggctacctgtatccagaaggagagaaaagaacagaggctggaccattggatcgaaggagacaaactccactggaaatgttc
tgtgaaggtaagtggtgcaaacagacattccagagagcctga