

PRR19-proline rich 19 Gene

Reference: PTXBC063021

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	proline-rich protein 19; proline rich 19		
For research use only.			

GENE DESCRIPTION	proline rich 19
GENE ID	284338
ORIGIN SPECIES	Human
ACCESSIONS	BC063021
VECTOR	pDONR223
NCBI SYMBOL	PRR19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggataccaggccagctctccagcctttcagcagcctgagaacctggtcgtgctcgtcgtcgaagactaggcgggaacgtaaacaggccctggtgggcagccgcccattagccaccacgacccctctgtggccattcgggatccacctgtggtccctactgcctccaag
ctcgtgtgtcataaccaggccggcgtgagccgggagcaccgggtctctcaaccacaggtgaaatccctagatgttgcaaggctgcttagcagtgaggacctggtgcccaggcagccccacactccccccaagccctcccccaagcccaggcagggcccaggaaccagccccac
gggtccagggaacaaagagaaccaggtgcctggaggttcgggcccaggcccacccaggtccccagagttgtctggtcgtggggcagctgtggcagagctgaggtgagctgagttgcccacaggccctccccggagggaacctgattcaggatgccaggatgccatcgtgcacacct
gcaggcctgtcatgtgtgtgctgacctgacctggtgcttcgggctgcccagccacctgcccaggtgagcgtccctcctgcccctgtactcccccttcttctgctggtcgtctagctcagctagcctctatgcttcttctctgtgccaggggccaagcctgggtctctga