

COX8A-cytochrome c oxidase subunit 8A (ubiquitous) Gene

Reference: PTXBC063025

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COX; COX8; COX8-2; COX8L; VIII; VIII-L; cytochrome c oxidase subunit 8A, mitochondrial; cytochrome c oxidase polypeptide VIII-liver/heart; cytochrome c oxidase subunit 8-2; cytochrome c oxidase subunit 8A (ubiquitous); cytochrome c oxidase subunit VIIIA (		

For research use only.

GENE DESCRIPTION	cytochrome c oxidase subunit 8A (ubiquitous)
GENE ID	1351
ORIGIN SPECIES	Human
ACCESSIONS	BC063025
VECTOR	pDONR223
NCBI SYMBOL	COX8A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgtcctgacgccgctgctgctcggggcttgacaggtcggccggggcgtccagtgccgcgcgccaaagatccattcgttgcgcggagggggaagctgggcatggaattggcgttgggcttacctcctgcttgcaccttccctcctgccaggggctggatcctgtcac
acctggagacctacaggaggccagagtga