

CXCL9-chemokine (C-X-C motif) ligand 9 Gene

Reference: PTXBC063122

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMK; Humig; MIG; SCYB9; crg-10; C-X-C motif chemokine 9; chemokine (C-X-C motif) ligand 9; gamma-interferon-induced monokine; monokine induced by gamma interferon; monokine induced by interferon-gamma; small-inducible cytokine B9; C-X-C motif chemokine li		

For research use only.

GENE DESCRIPTION	chemokine (C-X-C motif) ligand 9
GENE ID	4283
ORIGIN SPECIES	Human
ACCESSIONS	BC063122
VECTOR	pDONR223
NCBI SYMBOL	CXCL9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaaaagtgtgtcttttctcttggcgcacatcttgcgtgttctgattggagtgcaaggaaacccagtagtgagaaagggtcgtgttctgcacagcaccaacaaaggagactatccactacaatcttgaagaccttaacaattgccccaaagcccttctgcgagaaaattgaa
atcattgctacactgaagaatggagttcaaacatgtctaaacccagattcagcagatgtgaaggaaactgattaaaaagtgaggagaacagggtcagccaaaagaaaaagcaaaagaatgggaaaaaacatcaaaaaagaaagtctgaaagtctgaaaatctcaacgttctcgtcaaaaga
agactacataa