

C7orf46-chromosome 7 open reading frame 46 Gene

Reference: PTXBC063130

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C7orf46; C7orf46; protein FAM221A; family with sequence similarity 221 member A		

For research use only.

GENE DESCRIPTION	chromosome 7 open reading frame 46
GENE ID	340277
ORIGIN SPECIES	Human
ACCESSIONS	BC063130
VECTOR	pDONR223
NCBI SYMBOL	C7orf46
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcgggtgacgttgccctcgcggcgccggcggtggacgagctggagcaccggagaattgttggtgagatgatggaggaaacttttactcctgaagaatatgaagaatacaaaagaaagtgttacctcgcgttacaaaacagattattgtgagctggcggtcacc
aacagggatggattgtaactgtggccagagacactgtgtttgtacacataggtataaacaataaaactgacttggaaagcattctcagcagtgccccattgatctgccccgcaagtgcctgctccagtgccagggcttaccttatgtcccttgaatgtagccagccattcg
ctgcaggtgcaaacgctttgctgatcagcagtgctgcccgtgctttacatgcaatcatttcaagtggttcaggattccatagctgcttctactgtgcttgggtcagcctgcatatgccatgacacagtagtggaactaaagcaagaagattggctcagaaaaaccagtgaggacag
gacattccttatgcagccatggaggattaaatggtttcagctcgtggcggaaggctacatgctgttagatgacagtggttggtgtacattcagttgaatttttagaatctcccattacggcagtagacagccattcctaaaagcatttcaagcatcatctagtcttctccagaacgttaac
agatgtaggtacaagtagtcaagtttcttattaaggagacctgaagaggatgatgtgctttcttgaagacgataccaggaaagatgaaatggaaaaggctgctaagtggaaaggaaaagctccattgccatcagctacaaaaccttcatga