

PRRG4-proline rich Gla (G-carboxyglutamic acid) 4 (transmembrane) Gene

Reference: PTXBC063393

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRGP4; TMG4; transmembrane gamma-carboxyglutamic acid protein 4; proline rich Gla (G-carboxyglutamic acid) 4 (transmembrane); proline rich Gla 4 (transmembrane) isoform 1; proline rich Gla 4 (transmembrane) isoform 2; proline rich Gla 4 (transmembrane) is		

For research use only.

GENE DESCRIPTION	proline rich Gla (G-carboxyglutamic acid) 4 (transmembrane)
GENE ID	79056
ORIGIN SPECIES	Human
ACCESSIONS	BC063393
VECTOR	pDONR223
NCBI SYMBOL	PRRG4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttacgcttctgttctactcagccaactgccacagttaccctgggggttctcattgcgcaagaggtccaaggcttctaagcatcggggagaagaagtgtttacatcaaaagaagaagcaaaacttttcatacatagacgccttctgtataatagattgatctggagctcttcaactccgg
caacctagaagagagtgcaatgaagaactttgcaattatgaggaagccagagagatttttgggatgaagataaaacgattgcatttggcaggaatattcagctaaaggaccaaccacaaaatcagatggcaacagagagaaaatagattgtatggccttctgactggattaattgctgct
ggagtattttgttattttggattacttggctactatcttgtatcactaagtgtataggctacaacatccatgctcttcagccgtctatgaaagggggagggcacactccctccatcatttcagaagacctgaggaggctgccttgtctcattgccgcttctgtggaggatgcaggattaccttc
ttatgaacaggcagtggtggctgaccagaaaacacagtggttaccaccaccaccatctcgtggcacacaaaaggatttaggggtatttaaaaatctatgtctctccatctcactga