

SNRPF-small nuclear ribonucleoprotein polypeptide F Gene

Reference: PTXBC063397

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SMF; Sm-F; snRNP-F; small nuclear ribonucleoprotein F; sm protein F; small nuclear ribonucleoprotein polypeptide F		

For research use only.

GENE DESCRIPTION	small nuclear ribonucleoprotein polypeptide F
GENE ID	6636
ORIGIN SPECIES	Human
ACCESSIONS	BC063397
VECTOR	pDONR223
NCBI SYMBOL	SNRPF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtttaccctcaatcccaaccttctcaatggactaacaggaaagccagtgatggtgaaacttaagtggggaatggagtacaagggtatctggtatctgtagatggctacatgaacatgcagcttgcaaatacagaagaatacatagatggagcttgtctggacatctgggtgaagt
tttaataagggtgaataatgtcctttatcatcagaggtgtggaagaagaggagaagaatggggaatgagagaatag