

C21orf128-chromosome 21 open reading frame 128 Gene

Reference: PTXBC063412

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 21 open reading frame 128		
For research use only.			

GENE DESCRIPTION	chromosome 21 open reading frame 128
GENE ID	150147
ORIGIN SPECIES	Human
ACCESSIONS	BC063412
VECTOR	pDONR223
NCBI SYMBOL	C21orf128
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctggggctgcttgcacatcagaacactgctgtgccaacccacacctctcctgggtgctactcaacctccagcctacaggggctggagtatggggccagcgtggagatgccatgggaaacccggggtcctgcacgggtgagcttgagcctcagaccacacttccgcct
ggagagacacgatctgcattctcagcttcccacttgggtacaagtcagacaccactgtgtgggaaggagctctcagcctagccaagaagaggcaacagcagacatccatcaacttctcagaccattaagcagggttccaggtgtgaccgatggatggtgcttgagccatttctcctcta
caatcaagaggaggacacctgacgaccgacctctcagagcacgcagagaagtccgttcacagcaccttagctgggcttccctgggacggccggcctggcctggtgtgtgctggggactcacagtga