



Reference: PTXBC063435

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CFF; HCF; HCF-1; HCF1; HFC1; MRX3; PPP1R89; VCAF; host cell factor 1; VP16-accessory protein; protein phosphatase 1, regulatory subunit 89; host cell factor C1		

For research use only.

GENE DESCRIPTION	host cell factor C1 (VP16-accessory protein)
GENE ID	3054
ORIGIN SPECIES	Human
ACCESSIONS	BC063435
VECTOR	pDONR223
NCBI SYMBOL	HCFC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]