

RELL2-RELT-like 2 Gene

Reference: PTXBC063469

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C5orf16; RELT-like protein 2; receptor expressed in lymphoid tissues like 2; RELT like 2		

For research use only.

GENE DESCRIPTION	RELT-like 2
GENE ID	285613
ORIGIN SPECIES	Human
ACCESSIONS	BC063469
VECTOR	pDONR223
NCBI SYMBOL	RELL2
LENGTH	Proteogenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgtcgaaccacagcctgacctggaaccgccccacatggcctatatatgctctctctgtgtgtgtctctctctctatggcctgtgagctctatgatctgccagtgctcaagaagaagggtaccgtgccgcagctcaggggctctgagcctgacgatgccagctcagcccctgaggacgatgacatgaatgaggacacagtagagagatgttcgctgcatcaccagaatgaagccaatgctgagccttgaaggagatctggggacagtgaaaggagaaggacagctgacgtgtccagtgatgccactccagcctgacggacggagccccctccatcatcacacagctgacctgggcccctgacgcccctgcatccattgacccgacgcaagaggcctccactgtccgtcagggagcctccaaggaaggaaaaagccgccccggacaggggagaccactgtgtctctgtggcaggttccgggtgacacacattgagaagcgtatggactgcacgaacaccgtgatggctccccacagacagagctgggctctgtggggacagaccaggggtgtcaggggtctggggagggcagcccaaggcagggatgctgcatggagaggtgccccctgagaggccacagccccaggtcctagccagccccagtcagaatggaggactcagggacagcctaacccctgctgacttgaagggaacccagagctctgacagagccaactgaggccggaggaggcccaagcccaggctgccactcaagagggcaaatgggagccaagcaaacagacactctgatcaccaggtgctctaccacagggagcagggagtatgta
```