

SLC10A7-solute carrier family 10 (sodium/bile acid cotransporter family), member 7 Gene

Reference: PTXBC063471

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	sodium/bile acid cotransporter 7; Na(+)/bile acid cotransporter 7; SBF-domain containing protein; solute carrier family 10 (sodium/bile acid cotransporter family), member 7; solute carrier family 10 member 7		

For research use only.

GENE DESCRIPTION	solute carrier family 10 (sodium/bile acid cotransporter family), member 7
GENE ID	84068
ORIGIN SPECIES	Human
ACCESSIONS	BC063471
VECTOR	pDONR223
NCBI SYMBOL	SLC10A7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggctgctggagagaatgaggaagactgggtcatggtcgaatagtgctggcgcgctggagcctaaactggagccgtccataggggtgaatgggggaccactgaagccagaaataactgtatcctacattgctgttgaacaattctttaacagtggactatcattgaaaacaga
ggagctgaccagtgctttgggtcatctaaactgcattttttatcagatctttactcttgcattcttccagcaacaatatggctttttcttcagctttatcaatcacacccatcaacgaatggcttttaaagggttgcagacagtaggttgcattgcctccgctgtgtcttctgcagtgatttaacca
aggcagttgggtgaaatgaggcagctgcaatatattaatcagcctttggaagtgttttgaagtaaacatagtttaactgtctattacaacttttgcgtgga