

TM6SF1-transmembrane 6 superfamily member 1 Gene

Reference: PTXBC063487

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane 6 superfamily member 1		
For research use only.			

GENE DESCRIPTION	transmembrane 6 superfamily member 1
GENE ID	53346
ORIGIN SPECIES	Human
ACCESSIONS	BC063487
VECTOR	pDONR223
NCBI SYMBOL	TM6SF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgctctgcggccacccgggtctcgtgctgtccctcctggccatcccggtcacctatgtctcaaccacctggcgccagcatgattctggactattgtaggggtgtgctccctcatctgttctctgtagcactgctggctcgtgtcctcgtcaaaagaaaaccacccgggacc
cactgttctatgtatgcagttttggattaccagcgtggfgaacctcataggaactggagcaagatggaatcattgacgggttcacacactactgagagaggggaaccgtatctgaacaccgcataatgggcacatgatctgctactgggatgctcctcattatctgatgtacct
ggfcatggtggcagccatagcatgggaggaactatagaaccattggcctatattgggttgatctattatgagtggtgtgtttgtgccaggaacattgtagatgctggcatatattctattctgttctctactttgtgactgcactgtatggttagtggttctggtgatgtcctggtgcc
tga